




mixta

AR-EN



PATIENT TRANSFER TROLLEY

MHTS 9001

- Stretcher can be used in the intervention, hospital, and patient transfer.
- Back movement is made with the gas piston,
- Railings are polyurethane coated,
- Body structure is completely 304 quality stainless steel,
- Lying surface 8 mm thick compact laminate (HPL),
- 150 mm diameter wheels with brake system,
- IV pole is adjustable in height and has two hooks,
- There are 4 protective bumpers on each lower stretcher.



mixta

www.mixta.com.tr