
Genel Hastane Paslanmaz Çelik Ekipmanları / General Stainless Steel Hospital Equipments

TR / EN



MHTS 9000

- Stretcher can be used in emergency intervention, hospital and patient transfer.
- \bullet Stretcher can do two hydraulic hydraulic lifting movements and + / trendelenburg movements.
- The pedals on the lower right and left of the stretcher can make high-heeled movement and \pm / trendelenburg movements.
- Stretcher can make 80 degrees of backrest movement. The position can be achieved with a shock absorber under the pushing arm.
- The pushing handles of the stretcher are folded to prevent damage.
- The sides of the stretcher can be folded in one motion.
- The lying on the stretcher is made of compact laminate material.
- Stretcher has full size X-RAY cassette slot.
- Stretcher wheels 200 mm in diameter are locked in two and three ways.
- The mattress of the stretcher is made of an air permeable waterproof PU cover.
- The stretcher has four-cornered protective bumpers.
- All metal parts of the stretcher are painted with electrostatic powder paint.
- The lower trolley of the stretcher is covered with plastic material to store the oxygen tube and clothing.
- Dimensions: 2120x770x580/900Hmm





www.mixta.com.tr